ARTEFIT

SPCD0001

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TL.SPC.015 / Rev:12.03.2026/06

TABLE OF CONTENTS

1. Introduction
2. Intended Use
3. Warnings & Precautions
4. Intended Users
5. Hospital And Clinical Use
6. Home Use
7. Intended Patient Population
8. Contraindications
9. Device Description And Operating Principle
10. Instructions For Use
11. Therapy Profile Settings
12. Attaching And Donning The Garment
13. Operating The Device
14. Ending The Therapy
15. Technical Specifications
16. Garments And Accessories
17. Cleaning And Maintenance
18. Recommended Garments And Accessories For Use With The Product
19. Troubleshooting
20. Error Messages / Indicators
21. Warranty And Support
22. Device Service Life
23. Operating, Transport, And Storage Conditions
24. Hardware Specifications
25. Bluetooth Network Characteristics
26. Electromagnetic Compatibility
27. Conditions Requiring Consultation With Healthcare Professionals
28. Expected Clinical Benefit
29. Reporting Of Serious Incidents (Vigilance)
30. Device Label Symbol Glossary
31. Package Label Symbol Glossary

1- Introduction

Thank you for choosing the ARTEFIT SPCD0001 Sequential Pneumatic Compression Device, a professional-grade pneumatic compression system designed for both clinical and home use. Featuring an intuitive LCD touchscreen and compatible with both upper and lower extremity garments, the ARTEFIT SPCD0001 Sequential Pneumatic Compression Device delivers sequential gradient pressure therapy to support patients suffering from venous or lymphatic disorders.

2- Intended Use

The ARTEFIT SPCD0001 Sequential Pneumatic Compression Device is intended for the treatment and/or management of the following conditions:

- Lymphedema
- Chronic venous insufficiency
- Peripheral edema
- Postoperative edema
- Venous stasis ulcers
- As well as serving as a prophylactic device for Deep Vein Thrombosis (DVT).

This device is suitable for use in both clinical settings and home environments under medical supervision.

3 - Warnings & Precautions

- The information provided in this user manual is intended for general guidance on the use of the compression device and may not be applicable to every individual's specific clinical condition. Pressure levels, treatment duration, and frequency of application must be determined solely by the authorized healthcare professional(s) monitoring the patient.
- Pressure settings should not be changed unless explicitly recommended by the authorized healthcare professional. Application of high pressure levels should be evaluated with caution, as it may increase the risk of tissue damage.
- Under the following conditions, the device should be used only after careful evaluation and under the supervision of a healthcare professional:
 - Sensitive, irritated, sunburnt, or bruised skin, and skin with compromised integrity;
 - Skin conditions such as dermatitis, eczema, psoriasis, or skin cancer on or around the treatment area;
- If unusual changes such as blisters, redness, bumps, or discoloration are observed on the skin, or if sensations of burning, itching, or increased swelling are felt during or after treatment, the use of the device must be discontinued immediately and a healthcare professional must be consulted.
- Do not stand or walk while wearing the garments used with this device. Otherwise, there is a risk of injury due to tripping, slipping, or falling.
- This device is not intended for use during sleep. Use during sleep is not recommended, as it may adversely affect user safety.
- Do not carry or hang the device by its tubing, valve, or power cord.
- Keep the device away from CT or MR environments.
- Portable RF communications equipment (including peripherals such as mobile phones, antenna cables, and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device, including the power cord specified by the manufacturer. Otherwise, this could result in degradation of the equipment's performance.
- To avoid affecting the service life and normal function of the equipment, do not expose the device to direct sunlight and protect it from dirt and moisture.

- Operating the device near explosive or flammable gases may cause a fire hazard. Keep away from explosive, flammable, and volatile gases.
- Do not use in the presence of any contraindications.
- Avoid use during sleep.
- Do not use near water, in wet environments, or where aerosols are being sprayed.
- Do not operate the device on a soft surface, under a blanket or cover, or near a heat source. Inspect the device, accessories, and power cord for damage before use or cleaning.
- The garments are for single patient use only and must not be used by more than one person.
- For hygiene purposes and to prevent irritation, it is recommended to always wear clothing, bandages, or socks under the garments.
- Lie down in a comfortable position during operation. In addition, do not walk or stand while wearing the garments, as this may result in falls.
- Avoid folding, pinching, or kinking the tubing, as this may obstruct air flow.
- Do not wrap the tubing around your extremities, as this may restrict blood flow.
- Never attempt to open, repair, or modify the device; no modification of this equipment is allowed.
- Opening, repairing, or modifying the device by anyone other than authorized service personnel may adversely affect safety, performance, and regulatory compliance, and may pose a risk of injury.
- If you experience burning, itching, increased pain, numbness, or tingling, or observe changes in skin appearance such as discoloration, blisters, bumps, or increased swelling, discontinue use and consult a physician.
- Keep out of reach of children and pets.
- In the event of a power outage during treatment while the garments are inflated, first disconnect the garments from the device and allow air to escape from the garment connector. Then, unzip or open the touch-and-close fasteners and remove the garments from your body.
- To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective earth.
- This device must only be used with compatible compression garments, sold and provided separately by our company.
- There are no known side effects of the product.

4. Intended Users

- The device is designed for use in both home environments and hospital/clinical settings.
- In hospital and clinical settings, the setup, operation, parameter settings, and supervision of the device are performed by qualified healthcare providers (physicians, registered nurses, physiotherapists, and clinical technicians) in accordance with the instructions in the user manual; operational control of the device in these settings rests with healthcare professionals. Qualified healthcare providers, including physicians, registered nurses, physiotherapists, and clinical technicians, are responsible for the prescription, setup, operation, and supervision of the device in clinical settings.
- For home use, the setup and operation of the device are performed by the patient or caregiver, in accordance with the instructions in the user manual and the parameters recommended by qualified healthcare providers (physicians and physiotherapists). For home use, the patient or caregiver is responsible for operating the device, supervising its use, and entering the parameters provided by healthcare professionals.

- The device is also intended for use by patients with venous or lymphatic disorders. The device is suitable for use under medical supervision in the homes of patients diagnosed with Lymphedema, Chronic Venous Insufficiency, Peripheral Edema, Postoperative Swelling, Venous Stasis Ulcers, and Deep Vein Thrombosis (DVT) (prophylactically).
- Postoperative edema
- Venous stasis ulcers
- As well as serving as a prophylactic device for Deep Vein Thrombosis (DVT).
- The device is indicated for use both in clinical settings and in the homes of patients under medical supervision.

5. Hospital And Clinical Use

- The setup is performed by qualified healthcare providers (physicians, registered nurses, physiotherapists, and clinical technicians) in accordance with the instructions specified in the user manual.
- Device maintenance must be performed only by authorized service.
- Device Calibration: it is recommended to have the device calibrated once a year. Device calibration must be performed only by an authorized service provider.

6. Home Use

- The setup is performed by the patient or caregiver in accordance with the instructions in the user manual and the parameters recommended by qualified healthcare providers (physicians and physiotherapists).
- Device maintenance must be performed only by authorized service.
- Device Calibration: it is recommended to have the device calibrated once a year. Device calibration must be performed only by an authorized service provider.

7. Intended Patient Population

- The ARTEFIT Sequential Pneumatic Compression Device is intended for use with patients who have venous or lymphatic disorders, such as the following:
 - Lymphedema
 - Chronic venous insufficiency
 - Peripheral edema

- Based on FDA and CE standards, the ARTEFIT SPCD0001 Sequential Pneumatic Compression Device is suitable for use under medical supervision in patients within the 21-85 age group. There are no gender or weight restrictions within the patient population.

8. Contraindications

- Use is not recommended for individuals with the following conditions:
 - Active Deep Vein Thrombosis (DVT)
 - Suspected pulmonary embolism or clotting disorder
 - Limb infections (e.g., untreated cellulitis)
 - Congestive heart failure (unless medically controlled)
 - Lymphangiosarcoma
 - Untreated wounds, ulcers, or skin grafts
 - Severe peripheral artery disease
 - If you are pregnant, have impaired skin sensitivity, or have other vascular conditions, consult your physician before use.
 - Use outside the specified age groups is not recommended; consult your physician for exceptional cases.

9. Device Description And Operating Principle

- The Artefit Sequential Pneumatic Compression Device provides sequential, gradient pneumatic compression for the treatment of lymphedema and associated venous disorders. Sequential, gradient compression helps increase blood flow and facilitates the elimination of excess lymph by clearing it from the affected area of the body. The device provides sequential inflation/deflation cycles (from distal to proximal) at designated pressure levels.

• Adjustable / Monitorable Zone Pressure

- The pump pressure can be set between 20 and 120 mmHg. Individual pressure settings are possible for each zone prior to treatment, and real-time zone pressures can be monitored during the procedure.

• Monitorable Therapy Duration

- Therapy duration can be set to the following options: 15, 30, 45, and 60 minutes.
- The remaining time can be monitored during the therapy.

• Terminating the Therapy

- The therapy can be terminated at any time by selecting the 'Stop' option on the touchscreen.

• Pre-Therapy Pressure Application

- When using SPCG1 8-Chamber Proabdominal or SPCG8 16-Chamber Pro-pants garments, pre-therapy pressure may first be applied to the abdomen region under the supervision of a clinical technician, physiotherapist, or physician. The device then continues to operate at the set values.



• Extra Evacuation

- If extra evacuation of the air in the garment is required, the 'Extra Evacuation' option can be selected.

• Profile Memory

- The ARTEFIT Sequential Pneumatic Compression Device allows for the selection of 3 separate profiles. It can be adjusted for three different individuals, or three different programs can be entered for the same person. Each time a profile is selected, the last set values are automatically assigned.

• Front and Back Panel

1. LCD touchscreen
2. Garment Connection Port
3. Pin hole for aligning the garment socket
4. Power Cord Connection Socket
5. ON/OFF Button


10. Instructions For Use

• Preparing the Device for Use and Setup


- Remove the garment(s) from the box.
- Unfold the garment(s) and lay them out flat.
- Remove the device and the air mattress from the box – keep the packaging for transport and storage.
- Place the device on a flat and stable surface – the device should be close enough to be accessible to the controls during use.
- Ensure that the ambient temperature and humidity conditions are within the ranges defined in this manual.

• Connecting The Power Cord


- Insert the end of the power cord into the socket inlet on the back of the device.

- Plug the power cord into a grounded outlet.


- The device is ready for use; set the ON/OFF button to the ON (I) position.


• Configuring Device Settings on the LCD Touchscreen

- Select the device language by tapping the Language field on the screen.



• Selecting a Profile

- Select your username and/or program name from the LCD touchscreen.


- Save your username to one of the 3 profiles using the keyboard on the touchscreen. (Example profile: ABC)


• Selecting the Garment

- Once a profile is selected, the garment selection screen appears automatically.
- Select the product group recommended by clinical staff.


- Lower Body Group
 - 8 Chamber Proabdominal
 - 8 Chamber Leg Garment
 - 16 Chamber Pro-pants
- Upper Body Group
 - 8 Chamber Elite Pro Vest Garments

• Lower Body Group Garment Selection


- The screen displays default settings starting with 120 mmHg for Zone 1 and 85 mmHg for Zone 8. Select 'SETTINGS' on the touchscreen to adjust the values according to the instructions provided by clinical staff.

- If you are using the 8-Chamber Pro-abdominal or 16-Chamber Pro-pants garments, and the clinical staff recommends initial pre-therapy for the abdominal area, select 'ABDOMINAL ZONE' and set the desired duration. Once the treatment begins, the abdominal zone (Zones 7 and 8) operates for the set duration; then, the therapy starts from the feet towards the abdomen at the pre-adjusted pressure levels.

• Therapy Settings

- Press 'THERAPY' settings on the touch-screen.




- On the screen that appears:


- You can change the garment type.
- You can select the therapy duration as 15, 30, 45, or 60 minutes, according to the time recommended by clinical staff. The selected duration appears in green on the screen.
- You can save all selected settings into 3 different programs via Profile Settings.

• Default Settings (Presets 1, 2, 3)

- When an automatic setting is selected, the Default Settings are automatically configured as follows.
- Default 1 (Preset 1): Starts at 120 mmHg in Zone 1 (with a 5 mmHg gradual pressure decrease across zones).
- Default 2 (Preset 2): Starts at 75 mmHg in Zone 1 (with a 5 mmHg gradual pressure decrease across zones).
- Default 3 (Preset 3): Starts at 77 mmHg in Zone 1 (with a 3 mmHg gradual pressure decrease across zones).


- Select the recommended preset.
- For example, a 30-minute duration and Preset 1 have been selected below.

- Select 'Start' on the LCD touch-screen to begin the therapy.
- You can readjust your default settings and duration before starting by pressing the 'SETTINGS' tab.
- Profile settings change dynamically. The device memorizes the user's last program run and automatically starts with these settings for subsequent operations.


Upper Body Group Garment Selection

- Select the upper body group.
- After selection, all procedures continue as described for the lower body group.
- In the therapy settings for the upper body group, the Pre-Therapy Pressure Application is set as default and cannot be changed.


11. Therapy Profile Settings


- Select the profile you wish to run on the home page (Example: previously created ABC profile).

- Select 'SETTINGS' on the page that appears.


- Select GENERAL SETTINGS.


- If pre-therapy for the abdominal area is requested, select 'ABDOMINAL ZONE' and set the desired operation time. If selected, the ABDOMINAL ZONE button turns green.
- If abdominal pre-therapy pressure application is not requested, continue with the Therapy settings.


- Select AUTO SETTINGS.

- On the AUTO SETTINGS page, set the Zone 1 pressure value to the recommended level, and adjust the pressure change interval (DIFFERENTIAL VALUE) to the recommended setting. Press the SAVE button.


- Press 'START' to begin use. Your profile settings have been saved based on the values you entered. The next time you select your profile, you can directly start using it automatically with your last settings.


• Manual Settings

- Used to enter specific pressure values for each zone, from Zone 1 to Zone 8.
- Select the zone for which you wish to enter a value on the touchscreen.
- In the Differential Value section, select the increment interval for the pressure value of the chosen zone.
- In the Pressure Setting section, use the +/- buttons to adjust the pressure to your desired value.
- Press the SAVE button.
- The pressure value of each zone cannot be lower than the pressure value of the preceding zone. The system does not allow such an entry. Therefore, it is recommended to start from ZONE 1 and proceed upwards in sequence.


12. Attaching And Donning The Garment

- The garments are delivered ready for use, with the air tubing and the garment socket for device connection pre-mounted. It easily connects to the Garment Connection Port on the front panel of the device.
- While holding the garment socket in front of the device, ensure that the arrows on the socket are on the upper surface.
- While the arrows are on the upper surface, the pin of the garment socket must be on the right side. Ensure that the pin fits into the pin hole (3) on the garment connection port.

1. LCD Touchscreen
2. Garment Connection Port

Garment Socket Pin


3. Pin hole for aligning the garment socket

- Insert the tips and the pin of the garment socket into the holes and the pin hole of the garment connection port on the device.
- When using two garments, use the Air Distribution Manifold (ADM001). Following the arrows on the upper surface, insert the tips and the pin into the holes on the device. Then, insert the tips of both garment sockets into the holes of the Air Distribution Manifold (ADM001).

- Unzip or open the touch and-close fasteners of the garment with the tubing.
- Put on the garments, ensuring that the tubing runs along the sides of your legs or arms.


- Put on the Lower Body Garment so that the sole of your foot fits perfectly into the foot section of the garment.

- Put on the Upper Body Garment so that your fingertips reach the end of Zone 1.
- After putting on the garment, zip it up or close the touch-and-close fasteners.
- Do not kink the tubing.



• Using the AIR DISTRIBUTION MANIFOLD (ADM001)

- Based on the values recommended by qualified healthcare providers (doctors and physiotherapists), SPCG lower body garments can be used individually as right/left legs, or both right and left leg garments can be used simultaneously if recommended. In this case, the Air Distribution Manifold (ADM001) is used.

• **ADM001 Air Distribution Manifold Garment Connection Port**

- Insert the Right/Left Garment Sockets into the ADM001 connection port. Ensure that the arrow marks on the garment socket part of the ADM001 air distribution manifold are facing upwards, as shown in the figure.


- **ADM001 Air Distribution Manifold Garment Socket**
- **ADM001 Air Distribution Manifold Garment Socket Pin**

- The connection of the air distribution manifold to the SPCD001 device is performed as described in Article 12. Ensure that the arrow marks on the socket are facing upwards, make sure the garment socket pin fits into the socket pin hole (3) on the device, and insert the socket into the device port.
- Run the recommended program.
- The ADM001 air distribution manifold can only be used for the right and left 8-chamber leg garments in the lower body group.

13. Operating The Device


- Lean back comfortably within easy reach of the controls.
- Press the “On/Off” button to turn it on.
- Select your profile to access the home screen.
- Press the “Start” button on the screen and the device will begin to operate.
- The remaining therapy time can be tracked on the screen.
- The instantaneous pressure of the zones can be tracked individually on the screen. The current zone shown in the figure on the LCD touchscreen will be highlighted from light green to dark green as the desired pressure is reached.
- At the end of a timed therapy, the air inside the garment is deflated for 45 seconds. Once the evacuation process is complete, additional evacuation can be performed by answering “Yes” to the question “Should evacuation continue?” that appears on the screen. If the air inside the garment has been sufficiently deflated, select the “No” option. If the garment is to be folded and packed after therapy, additional evacuation is recommended.

- Once the process is complete, turn off the device using the “On/Off” button.
- To disconnect the garment from the device, pull and remove the garment socket.
- Unzip or open the touch-and-close fasteners and remove the garment.

14. Ending The Therapy

- Select the “Stop” option on the LCD touchscreen to end the therapy.
- The evacuation process will start immediately, and air evacuation will continue for 45 seconds.

- Select either the “Yes” or “No” option for additional evacuation, depending on the need to deflate the garment.
- If “No” is selected, the screen returns to the home page.

15. Technical Specifications

TITLE	VALUE
Product Name	SEQUENTIAL PNEUMATIC COMPRESSION DEVICE
Power Supply Rating	110–220 V AC, 40/50 Hz
Nominal Input	12 V DC, 4 A
Electrical Classification	Class II
Applied Part	Type B
Ingress Protection	IP21
Power Supply Isolation	By unplugging
Mode of Operation	Continuous
Essential Performance	Cyclical inflation and deflation of the garment(s) by the device
Treatment Duration	15 / 30 / 45 / 60 minutes
Pressure Range	20 – 120 mmHg
Precision	3 mmHg
Accuracy	± %20
Functions	• Adjustable / Monitorable Zone Pressure • Monitorable Therapy Duration • Terminating the Therapy • Pre-Therapy Pressure Application • Extra Evacuation • Profile Memory

Warranty	Device: 3 years / Garment: 1 year
Expected Service Life	3 years
Software Safety Class	A
Regulatory Classification	EU MDR 2017/745 – Class IIa
Weight	2,4 kg
Dimensions (HxWxL)	117,3 × 268,3 × 245 mm

16. Garments And Accessories

The ARTEFIT SPCD0001 Sequential Pneumatic Compression Device is fully and exclusively compatible with Artefit garments. Artefit garments are designed for targeted and comfortable therapy and are available for both upper and lower extremities.

• Garment Types:

- Upper Body Group: Designed to cover the entire arm. Available in various sizes and chamber configurations.
- Lower Body Group: Covers from the foot up to the upper thigh. Designed to provide compression along the limb.
- Bilateral Use: In the lower body group, single-leg garments can be used simultaneously for both the right and left sides. The air distribution manifold supports dual-garment therapy.

• Key Features:

- Multi-chamber design for sequential gradient inflation
- Durable and easy-to-clean material
- Zippers and/or touch-and-close fasteners for a secure and adjustable fit

• Included Accessories:

- Power Adapter: Medical-grade and region-specific plug

17. Cleaning And Maintenance

Proper maintenance extends the life of the ARTEFIT Sequential Pneumatic Compression Device and ensures hygienic therapy.

• Daily Maintenance:

- Wipe the exterior of the device with a soft, damp cloth.
- Do not use abrasive cleaners or solvents, and do not immerse the device in water.

• Garment Cleaning:

- Disconnect the garment from the device.
- Unfasten the garment completely and lay it flat on a level surface.
- Gently wipe with a damp and/or soapy cloth.
- Do not machine wash, do not tumble dry, or do not immerse in water.
- Wait for the garment to dry completely before the next use.

• Tubing Maintenance:

- Ensure there are no kinks or blockages in the tubing.
- Wipe with a dry cloth if necessary.
- Do not bend the tubing sharply.

• Storage:

- Store in a cool, dry place away from direct sunlight.
- Avoid placing heavy objects on the device or garments.
- Keep away from moisture, heat sources, or open flames.

- **Service:**

- There are no user-serviceable parts inside the device.
- Contact your authorized service provider for service, repairs, or technical issues.
- Perform regular checks for signs of wear or damage. If you detect any abnormality, discontinue use and seek support.
- At the end of their service life, the device and its accessories must be disposed of in accordance with applicable national legislation and environmental regulations.

18. Recommended Garments And Accessories For Use With The Product

The ARTEFIT SPCD0001 Sequential Pneumatic Compression Device must only be used with compatible SPCG compression garments, which are provided and sold separately by our company.

There are no known restrictions on combinations of devices and equipment.

The use of garments manufactured from biocompatible materials is recommended with the ARTEFIT SPCD0001 Sequential Pneumatic Compression Device.

- **Garments**

MODEL/REF NO	PRODUCT NAME	BODY SIZE
SPCD1-SPCG10001	8 chamber Proabdominal garments	Small L
SPCD1-SPCG10002	8 chamber Proabdominal garments	Small R
SPCD1-SPCG10003	8 chamber Proabdominal garments	Medium L
SPCD1-SPCG10004	8 chamber Proabdominal garments	Medium R
SPCD1-SPCG10005	8 chamber Proabdominal garments	Large L
SPCD1-SPCG10006	8 chamber Proabdominal garments	Large R
SPCD1-SPCG40001	8 chamber leg garments	Small S
SPCD1-SPCG40002	8 chamber leg garments	Medium S
SPCD1-SPCG40003	8 chamber leg garments	Large S
SPCD1-SPCG40004	8 chamber leg garments	Small N
SPCD1-SPCG40005	8 chamber leg garments	Medium N
SPCD1-SPCG40006	8 chamber leg garments	Large N
SPCD1-SPCG40007	8 chamber leg garments	Small W
SPCD1-SPCG40008	8 chamber leg garments	Medium W
SPCD1-SPCG40009	8 chamber leg garments	Large W
SPCD1-SPCG40010	8 chamber leg garments	Medium EW
SPCD1-SPCG60001	8 chamber Elite Pro Vest garments	Small L

SPCD1-SPCG60002	8 chamber Elite Pro Vest garments	Small R
SPCD1-SPCG60003	8 chamber Elite Pro Vest garments	Medium L
SPCD1-SPCG60004	8 chamber Elite Pro Vest garments	Medium R
SPCD1-SPCG60005	8 chamber Elite Pro Vest garments	Large L
SPCD1-SPCG60006	8 chamber Elite Pro Vest garments	Large R
SPCD1-SPCG70001	16 chamber elite bilateral garments	Small L
SPCD1-SPCG70002	16 chamber elite bilateral garments	Small R
SPCD1-SPCG70003	16 chamber elite bilateral garments	Medium L
SPCD1-SPCG70004	16 chamber elite bilateral garments	Medium R
SPCD1-SPCG70005	16 chamber elite bilateral garments	Large L
SPCD1-SPCG70006	16 chamber elite bilateral garments	Large R
SPCD1-SPCG80001	16 chamber pro pants	Small
SPCD1-SPCG80002	16 chamber pro pants	Medium
SPCD1-SPCG80003	16 chamber pro pants	Large

Abbreviations Table:

- L : Left
- R : Right
- S : Small S
- N: Narrow
- W: Wide
- EW : Extra Wide

• **Air Distribution Manifold**




MODEL NO/REF NO	PRODUCT NAME
SPCD1-ADM001	AIR DISTRIBUTION MANIFOLD

19. Troubleshooting

Issue	Potential Cause	Solution
Device does not turn on.	Power cord is not plugged in or the outlet is not active.	Check the power connection and verify that the outlet is working.
Touchscreen does not respond.	Software lag or display failure.	Turn the device off and on again; if the problem persists, contact the support team
No air pressure in the garment	Tubing connection is disconnected or blocked.	Ensure all connections are secure and the tubing is not kinked
Uneven pressure in chambers	Garment is not positioned correctly or there is a chamber leak.	Re-apply the garment and check for visible damage.
Excessive noise or vibration	Internal compressor issue	Discontinue use and contact an authorized service center.
Error message on screen	Internal system error	Discontinue use and contact support.

20. Error Messages / Indicators

The Artefit SPCD0001 Sequential Pneumatic Compression Device does not include an audible alarm system. Visual warning indicators are displayed on the LCD screen and can be referenced in the user manual.

- **E01** - High Pressure Error Indicator: This may occur if there is a kink or other blockage in the tubing connecting the garments, or if the connection port is blocked. The E01 message is displayed on the LCD touchscreen.
- **E02** - Low Pressure Error Indicator: This may occur if the garment is not properly secured to the connection port or if there is a hole in the garment causing a leak. The E02 message is displayed on the LCD touchscreen.
- **E03** - Software System Error Indicator: This may occur if an internal software system error occurs that causes the device to fail to operate. The E03 message is displayed on the LCD touchscreen.

• Additional Warnings:

- Ensure that the garment zippers/touch-and-close fasteners are completely closed/secured before each session.
- Avoid using extension cords or power strips.
- After 45 and 60-minute sessions, allow the device to rest for 5 minutes before the next session.
- If problems persist after following the steps above, discontinue use and contact the user support team listed in the next section.

21. Warranty And Support

The ARTEFIT SPCD0001 Sequential Pneumatic Compression Device comes with a limited warranty covering defects in materials and workmanship. Please keep your purchase receipt. For support, troubleshooting, or replacement parts, please contact the customer service center listed below. The ARTEFIT SPCD0001 Sequential Pneumatic Compression Device is maintenance-free during the warranty period, and annual calibration is recommended. Devices sent to the manufacturer for service/repair are returned to the customer after post-service performance and safety validations are performed in accordance with the technical service instructions.

Product Name / Ref. No. : Sequential Pneumatic Compression Device / SPCD0001

Manufacturer: Sanko Tekstil İşletmeleri
 San. Ve Tic. A.Ş. İSKO Şubesi
 Address: Organize Sanayi Bölgesi 3. Cad.
 İnegöl / Bursa / Türkiye
 Support Phone Number: +902242807878
 E-Mail: info@artefit.com

22. Device Service Life

The service life of the ARTEFIT SPCD0001 Sequential Pneumatic Compression Device is 3 years. The service life of ARTEFIT SPCG Sequential Pneumatic Compression Garments, compatible with the SPCD0001 device, is 1 year.

23. Operating, Transport, And Storage Conditions

- Emissions generated during normal use
 - Compressed air
 - Minimum acoustic energy, quiet operation
- Operating Environment
 - Designed for use in healthcare facilities or home environments.
 - Not suitable for use in CT and MRI environments, or in environments rich in explosive, flammable, and volatile gases.
 - Suitable for use in dry environments.
 - Temperature: 00C - 600C
 - Relative Humidity: 10% - 70% RH
 - Atmospheric Pressure: 700 - 1060 hPa

24. Hardware Specifications

- The device features a high-performance 32-bit microprocessor as its minimum hardware architecture and a touchscreen interface for user interaction.
- A pressure sensor, valve driver units, and pump motor control circuits are incorporated in the system as critical hardware components for the software to execute its core functions.
- The system operates stably on a mainboard powered by 12 Volts and incorporating the necessary power management units.

25. Bluetooth Network Characteristics

- The device is not intended to connect to any Bluetooth (BT) network (IT-network) or to transfer data via external networks. The device does not connect to any BT network or transfer data via external networks.
- Since the device has no physical connection ports accessible to the outside world, such as USB or RS485, the requirements regarding BT network characteristics are not applicable.
- Internal data communication is performed exclusively as a closed-circuit via a serial communication protocol between the device's main processor and the display controller.

26. Electromagnetic Compatibility

Emissions	Compatibility	Electromagnetic Environment - Guidance
RF emissions	Group 1	The device uses RF energy only for its internal functions. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions	POSITIVE	The device complies with the harmonic current limits specified in the standard.
Voltage fluctuations and flicker emissions	POSITIVE	The device complies with the voltage fluctuation and flicker limits.

• Electromagnetic Emissions

Immunity Test	Immunity Test Level
IEC 61000-4-2 Electrostatic Discharge (ESD) Immunity	±8 kV contact, ±2, 4, 8, 15 kV air discharge
IEC 61000-4-3 Radiated RF Electromagnetic Field Immunity	80 MHz – 2.7 GHz, 10 V/m, 80% AM at 1 kHz
IEC 61000-4-3 Proximity fields from RF wireless communications equipment -Immunity	TS EN 60601-1-2:2016/A1:2021 Standard Table 9 (9 V/m - 28 V/m)
IEC 61000-4-4 Electrical Fast Transient / Burst	±2 kV Power Lines (100 kHz)
IEC 61000-4-5 Surge Immunity	±0.5 kV – ±1 kV Line -to-Line
IEC 61000-4-6 Conducted RF Immunity	150 kHz – 80 MHz, 3 V – 6 V (rms), 80% AM at 1 kHz
IEC 61000-4-8 Power Frequency Magnetic Field Immunity	30 A/m, 50 Hz or 60 Hz
IEC 61000-4-11 Voltage Dips	0% Ut for 0.5 cycle, 0% Ut for 1.0 cycle, 70% Ut for 25/30 cycles
IEC 61000-4-11 Voltage Interruptions	0% Ut for 250/300 cycles

27. Conditions Requiring Consultation With Healthcare Professionals

If any of the following situations occur while using the device at home, use of the device should be discontinued immediately, and a healthcare professional should be consulted:

- Redness, irritation, sores, or changes in skin color
- Increased pain during or after use
- Numbness, tingling, or loss of sensation
- Increased swelling in the area of application
- Sudden shortness of breath, chest pain, or dizziness

28. Expected Clinical Benefit

This device is designed to support venous and lymphatic circulation in the lower and upper extremities by operating on the principle of sequential pneumatic compression. When used properly and consistently, it is expected to help reduce edema, alleviate feelings of fullness and tension in the extremities, and enhance user comfort.

ARTEFIT

Product Name: Sequential Pneumatic Compression Device

Ref.: SPCD0001 UDI-DI:08684523754033
Manufacturer and Authorized Service:
Sanko Tekstil İşletmeleri San. ve Tic. A.Ş.
ISKÖ Division

Address: Organize Sanayi Bölgesi 3. Cad.
İnegöl / Bursa / Türkiye
Phone: +902242807878
Email: info@artefit.com
Notified Body: Kiwa Belgelendirme
Hizmetleri A.Ş.

Address: İTOSB 9. Cadde No 15 Tepeören
Tuzla – İstanbul – Türkiye
Phone: +90 (216) 593 25 75
E-posta: posta@kiwa.com.tr

29. Reporting Of Serious Incidents (Vigilance)

The user, caregiver, and/or healthcare professional must report any serious adverse event (e.g., death, unexpected serious deterioration in health, hospitalization, permanent injury, life-threatening condition, etc.) occurring during the use of this device or related to it to the manufacturer without delay. Prompt reporting is critical for the evaluation of the event and the initiation of necessary corrective and preventive actions.

MANUFACTURER NOTIFICATION CHANNELS:

- Phone: +902242807878
- Email: info@artefit.com
- Address: Organize Sanayi Bölgesi 3. Cad. İnegöl / Bursa / Türkiye

Contents of the report (if possible): device/garment name/reference/model, serial/lot/UDI information, date and location of the event, description of the event, condition of the affected person, any photos/additional documents, and your contact information.

Serious adverse events: The event must be reported to us by the user, caregiver, healthcare professional, or affected person(s) within 10 days of its occurrence.

Death or unexpected serious deterioration in health: The event must be reported to us by the user, caregiver, healthcare professional, or affected person(s) within 7 days of its occurrence.

Serious public health threat: The event must be reported to us by the user, caregiver, healthcare professional, or affected person(s) within 1 day of its occurrence.

NOTIFICATION TO THE COMPETENT AUTHORITY:

Users, caregivers, and healthcare professionals may, if they wish, report a serious event directly to the competent authority in addition to notifying the manufacturer..

Competent Authority for Türkiye:

Türkiye

MDR-IVDR

Turkish Medicines and Medical Devices Agency (TITCK) Ministry of Health

Turkish Medicines and Medical Devices Agency

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